

解决方案 | 原油市场中的重金属检测

原油是一种粘稠的、深褐色（有时有点绿色的）液体，是石油刚开采出来未经提炼或加工的物质。主要由碳、氢、硫、氮、氧五种元素组成（质量分数）。此外，石油中还含有微量的钙、镁、铁、镍、铜、钒、砷、铝、硅等金属元素。这些金属元素虽然含量很少，但对炼油生产及催化剂失活等影响极大。

金属元素的检测方法主要有原子发射光谱法、原子吸收光谱法、高效液相色谱法，电化学分析法等。其中电感耦合等离子体发射光谱法，因其具有操作简便、快速且可同时测定多种元素的优点，受到分析工作者的青睐。

本文利用 GBC Quantima ，湿法酸消解对原油中的金属元素进行检测，检测结果表明，该方法便捷、高效、准确度高，可以满足原油、炼油厂及下游装置的质量控制需求。

实验步骤：称取 1.0g（精确到万分之一）样品于烧杯中，加入 25mL 硝酸和 4mL 高氯酸，加盖表面皿，电热板 150 度低温加热 30min，调整温度至 200 度继续加热至溶液呈淡黄色液体，取下冷却至室温，开盖，继续加热蒸至湿盐状，再加入 3mL 硝酸溶解，溶解后取下冷却，用去离子水转移定容至 25-50mL，摇匀，备用。

详细报告

Method Name : 3

El	nm	Type	Mono	Order	Power (W)	Neb Height (L/min) (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT (V)	Int. (s)	Reps	Bkgd Corr
Pb	II 220.353	An	1	2	1000	0.80 8.0	10.0	0.5	10	600	0.50	3	Fixed

Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1


Pb	Conc = 4.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 14413.5
	Net Int %RSD = 2.69	Peak Int = 15639.7		

Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Pb	Conc = 2.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 7016.24
	Net Int %RSD = 3.24	Peak Int = 8247.67		

Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1				
Pb	Conc = 1.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 3567.74
	Net Int %RSD = 5.08	Peak Int = 4744.33		
Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1				
Pb	Conc = 0.5000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 1719.34
	Net Int %RSD = 5.25	Peak Int = 2998.33		
Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1				
Pb	Conc = 0 mg/L	Conc %RSD =	Flags = b	Net Int = 28.2572
	Net Int %RSD = 376	Peak Int = 1454.33		

Pb II 220.353 nm (mg/L)



Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	14413.5	4.000	4.018	0.440
Standard 2	7016.24	2.000	1.963	-1.85
Standard 3	3567.74	1.000	1.005	0.526
Standard 4	1719.34	0.5000	0.4919	-1.62
Blank	28.2572	0	0.02219	***

Calibration Coefficients	C0	C1	R
	0.014343	2.777e-4	0.9998


Sample Blank [Sample Blank] Act Wt = 1 Nom Wt = 1 Dil = 1				
Pb	Conc = 0.02944 mg/L	Conc %RSD = 180	Flags = b	Net Int = 54.3430
	Net Int %RSD = 351	Peak Int = 1244.33		
1-2 Act Wt = 1.6292 Nom Wt = 25 Dil = 1				
Pb	Conc = -0.2410 mg/L	Conc %RSD = 671	Flags =	Net Int = -2.197
	Net Int %RSD = 189	Peak Int = 1119.00		
1-1 Act Wt = 1.455 Nom Wt = 25 Dil = 1				
Pb	Conc = -0.7931 mg/L	Conc %RSD = -76	Flags = b	Net Int = -111.846
	Net Int %RSD = -1e2	Peak Int = 1115.67		
2-1 Act Wt = 1.401 Nom Wt = 25 Dil = 1				

Pb	Conc = -0.5898 mg/L %RSD = -2e2	Conc %RSD = -1e2 Peak Int = 1069.67	Net Int	Flags = b	Net Int = -64.6688
2-2 Act Wt = 1.5466 Nom Wt = 25 Dil = 1					
Pb	Conc = -1.146 mg/L %RSD = -66	Conc %RSD = -52 Peak Int = 1041.67	Net Int	Flags =	Net Int = -200.965
3-1 Act Wt = 1.495 Nom Wt = 25 Dil = 1					
Pb	Conc = -0.6633 mg/L %RSD = -2e2	Conc %RSD = -1e2 Peak Int = 1162.33	Net Int	Flags = b	Net Int = -88.4759
3-2 Act Wt = 1.3944 Nom Wt = 25 Dil = 1					
Pb	Conc = -0.6444 mg/L %RSD = -32	Conc %RSD = -18 Peak Int = 1192.33	Net Int	Flags =	Net Int = -75.0675

Method Name :3

El	nm	Type	Mono	Order	Power (W)	Neb (L/min)	Height (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT (V)	Reps	Bkgd Corr
Na I	588.995	An	1	1	1000	0.80	6.0	10.0	0.5	10	500	3	Fixed
Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1													
Na	Conc = 5.000 mg/L		Conc %RSD = 0.00		Flags =		Net Int = 192108						
	Net Int %RSD = 1.11		Peak Int = 200828										
Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1													
Na	Conc = 2.000 mg/L		Conc %RSD = 0.00		Flags =		Net Int = 81421.0						
	Net Int %RSD = 0.20		Peak Int = 90178.0										
Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1													
Na	Conc = 1.000 mg/L		Conc %RSD = 0.00		Flags =		Net Int = 43484.4						
	Net Int %RSD = 1.20		Peak Int = 52272.7										
Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1													
Na	Conc = 0.5000 mg/L		Conc %RSD = 0.00		Flags =		Net Int = 24120.1						
	Net Int %RSD = 2.59		Peak Int = 32587.3										
Standard 5 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1													
Na	Conc = 0.1000 mg/L		Conc %RSD = 0.00		Flags =		Net Int = 8821.11						
	Net Int %RSD = 2.72		Peak Int = 17264.0										
Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1													
Na	Conc = 0 mg/L		Conc %RSD =		Flags = b		Net Int = 4168.95						
	Net Int %RSD = 3.18		Peak Int = 13468.7										

Na I 588.995 nm (mg/L)



Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	192108	5.000	4.983	-0.335
Standard 2	81421.0	2.000	2.031	1.54
Standard 3	43484.4	1.000	1.019	1.89
Standard 4	24120.1	0.5000	0.5024	0.482
Standard 5	8821.11	0.1000	0.09433	-5.67
Blank	4168.95	0	-0.02976	***

Calibration Coefficients	C0	C1	R
	-0.14096	2.667e-5	0.9999

Sample Blank [Sample Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

Na Conc = 0.6925 mg/L Conc %RSD = 0.68 Flags = b Net Int = 31246.0
Net Int %RSD = 0.57 Peak Int = 41204.7

1-1 Act Wt = 0.8594 Nom Wt = 50 Dil = 2

Na Conc = 51.94 mg/L Conc %RSD = 3.24 Flags = Net Int = 47980.1
Net Int %RSD = 3.77 Peak Int = 56578.0

1-2 Act Wt = 0.7595 Nom Wt = 50 Dil = 2

Na Conc = 51.41 mg/L Conc %RSD = 2.88 Flags = Net Int = 45884.3
Net Int %RSD = 1.87 Peak Int = 54825.3

2-1 Act Wt = 1.1685 Nom Wt = 50 Dil = 2

Na Conc = 21.21 mg/L Conc %RSD = 2.80 Flags = Net Int = 40535.9
Net Int %RSD = 1.10 Peak Int = 48992.0

2-2 Act Wt = 1.4356 Nom Wt = 50 Dil = 2

Na Conc = 21.09 mg/L Conc %RSD = 2.92 Flags = Net Int = 42594.9
Net Int %RSD = 2.82 Peak Int = 51894.0

3-1 Act Wt = 1.086 Nom Wt = 50 Dil = 2

Na Conc = 19.96 mg/L Conc %RSD = 2.14 Flags = Net Int = 39374.9
Net Int %RSD = 2.02 Peak Int = 47895.3

3-2 Act Wt = 1.2581 Nom Wt = 50 Dil = 2

Na Conc = 21.83 mg/L Conc %RSD = 3.00 Flags = Net Int = 41542.2
Net Int %RSD = 2.922 Peak Int = 50174.7

Method Name :3

El	nm	Type	Mono	Order	Power (W)	Neb (L/min)	Height (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT Int. (V)	Int. (s)	Reps	Bkgd Corr	
V	II	311.838	An	1	1	1000	0.80	6.0	10.0	0.5	10	650	0.50	3	Fixed

Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

V	Conc = 2.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 125422
	Net Int %RSD = 2.50	Peak Int = 133001		

Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

V	Conc = 1.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 65580.1
	Net Int %RSD = 2.93	Peak Int = 72892.7		

Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

V	Conc = 0.5000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 30547.6
	Net Int %RSD = 1.96	Peak Int = 37520.0		


Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

V	Conc = 0.1000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 6778.09
	Net Int %RSD = 4.95	Peak Int = 13830.7		

Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

V	Conc = 0 mg/L	Conc %RSD =	Flags = b	Net Int = -93.3638
	Net Int %RSD = -12	Peak Int = 6882.00		

V II 311.838 nm (mg/L)



Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	125422	2.000	1.986	-0.719
Standard 2	65580.1	1.000	1.036	3.62
Standard 3	30547.6	0.5000	0.4805	-3.91
Standard 4	6778.09	0.1000	0.1033	3.35
Blank	-93.3638	0	-0.005666	***

Calibration Coefficients	C0	C1	R
	-0.00418	1.586e-5	0.9993

Sample Blank [Sample Blank] Act Wt = 1 Nom Wt = 1 Dil = 1				
V	Conc == -0.008789 mg/L Net Int %RSD = -2e2	Conc %RSD = -1e2 Peak Int = 6759.33	Flags = b	Net Int = -290.217
1-1 Act Wt = 1.4546 Nom Wt = 25 Dil = 1				
V	Conc = 0.6808 mg/L Net Int %RSD = 7.63	Conc %RSD = 7.32 Peak Int = 8571.33	Flags=b	Net Int = 2206.61
1-2 Act Wt = 1.6292 Nom Wt = 25 Dil =				
V	Conc = 0.6360 mg/L Net Int %RSD = 10.1	Conc %RSD = 8.97 Peak Int = 8922.67	Flags = b	Net Int = 2322.11
2-1 Act Wt = 1.4015 Nom Wt = 25 Dil = 1				
V	Conc = 1.883 mg/L Net Int %RSD = 2.86	Conc %RSD = 2.73 Peak Int = 13048.0	Flags = b	Net Int = 6364.76
2-2 Act Wt = 1.5466 Nom Wt = 25 Dil = 1				
V	Conc = 1.838 mg/L Net Int %RSD = 2.25	Conc %RSD = 2.52 Peak Int = 13718.0	Flags =	Net Int = 6878.35
3-1 Act Wt = 1.495 Nom Wt = 25 Dil = 1				
V	Conc = 5.473 mg/L Net Int %RSD = 3.19	Conc %RSD = 3.14 Peak Int = 26582.7	Flags = b	Net Int = 20339.6
3-2 Act Wt = 1.3944 Nom Wt = 25 Dil = 1				
V	Conc = 5.389 mg/L Net Int %RSD = 3.27	Conc %RSD = 3.32 Peak Int = 25166.0	Flags =	Net Int = 18655.1

Method Name :3


El	nm	Type	Mono	Order	Power (W)	Neb (L/min)	Height (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT (V)	Int. (s)	Reps	Bkgd Corr
Al I	396.152	An	1	1	1000	0.80	6.0	10.0	0.5	10	500	0.50	3	Fixed
Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1														
Al	Conc = 8.000 mg/L Net Int %RSD = 0.57		Conc %RSD = 0.00 Peak Int = 250679		Flags =		Net Int = 237190							
Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1														
Al	Conc = 4.000 mg/L Net Int %RSD = 0.13		Conc %RSD = 0.00 Peak Int = 133219		Flags =		Net Int = 119890							
Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1														
Al	Conc = 2.000 mg/L Net Int %RSD = 2.25		Conc %RSD = 0.00 Peak Int = 72880.7		Flags =		Net Int = 59725.7							
Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1														
Al	Conc = 1.000 mg/L Net Int %RSD = 2.07		Conc %RSD = 0.00 Peak Int = 44546.7		Flags =		Net Int = 31654.9							
Standard 5 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1														

Al	Conc = 0.5000 mg/L Net Int %RSD = 2.34	Conc %RSD = 0.00 Peak Int = 27805.3	Flags =	Net Int = 14951.3
----	---	--	---------	-------------------

Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

Al	Conc = 0 mg/L Net Int %RSD = 14.9	Conc %RSD = Peak Int = 13115.3	Flags = b	Net Int = -393.601
----	--------------------------------------	-----------------------------------	-----------	--------------------

Al I 396.152 nm (mg/L)



Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	237190	8.000	7.987	-0.159
Standard 2	119890	4.000	4.023	0.569
Standard 3	59725.7	2.000	1.989	-0.534
Standard 4	31654.9	1.000	1.041	4.06
Standard 5	14951.3	0.5000	0.4760	-4.79
Blank	393.601	0	-0.01600	***

Calibration Coefficients	C0	C1	R
	-0.02930	3.380e-5	0.9999

Sample Blank [Sample Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

Al	Conc = -0.01329 mg/L Net Int %RSD = 9.75	Conc %RSD = 11.6 Peak Int = 13176.0	Flags = b	Net Int = 473.804
----	---	--	-----------	-------------------

1-1 Act Wt = 0.8595 Nom Wt = 50 Dil = 2

Al	Conc = 145.0 mg/L Net Int %RSD = 2.66	Conc %RSD = 2.71 Peak Int = 50290.7	Flags = b	Net Int = 37344.7
----	--	--	-----------	-------------------

1-2 Act Wt = 0.7595 Nom Wt = 50 Dil = 2

Al	Conc = 149.3 mg/L Net Int %RSD = 2.31	Conc %RSD = 2.36 Peak Int = 47451.3	Flags = b	Net Int = 34014.5
----	--	--	-----------	-------------------

2-1 Act Wt = 1.4356 Nom Wt = 50 Dil = 2

Al	Conc = 3.847 mg/L Net Int %RSD = 3.86	Conc %RSD = 3.84 Peak Int = 15006.7	Flags = b	Net Int = 2107.72
----	--	--	-----------	-------------------

2-2 Act Wt = 1.1685 Nom Wt = 50 Dil = 2

Al	Conc = 4.403 mg/L Net Int %RSD = 3.86	Conc %RSD = 3.68 Peak Int = 14417.3	Flags =	Net Int = 1996.22
----	--	--	---------	-------------------

3-1 Act Wt = 1.0860 Nom Wt = 50 Dil = 2

Al	Conc = 12.54 mg/L Net Int %RSD = 2.58	Conc %RSD = 2.79 Peak Int = 17316.0	Flags = b	Net Int = 4508.18
----	--	--	-----------	-------------------

3-2 Act Wt = 1.2581 Nom Wt = 50 Dil = 2

Al	Conc = 12.19 mg/L Net Int %RSD = 2.29	Conc %RSD = 2.89 Peak Int = 17440.0	Flags =	Net Int = 5015.62
----	--	--	---------	-------------------

Method Name :3

El	nm	Type	Mono	Order	Power (W)	Neb (L/min)	Height (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT (V)	Int. (s)	Reps	Bkgd Corr	
Ni	II	231.604	An	1	2	1000	0.80	6.0	10.0	0.5	10	550	0.50	3	Fixed

Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Ni	Conc = 4.000 mg/L Net Int %RSD = 0.57	Conc %RSD = 0.00 Peak Int = 85982.0	Flags =	Net Int = 83753.9
----	--	--	---------	-------------------

Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Ni	Conc = 2.000 mg/L Net Int %RSD = 2.81	Conc %RSD = 0.00 Peak Int = 45055.3	Flags =	Net Int = 43082.7
----	--	--	---------	-------------------

Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Ni	Conc = 1.000 mg/L Net Int %RSD = 2.36	Conc %RSD = 0.00 Peak Int = 24982.0	Flags =	Net Int = 23140.3
----	--	--	---------	-------------------


Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Ni	Conc = 0.5000 mg/L Net Int %RSD = 1.47	Conc %RSD = 0.00 Peak Int = 12684.7	Flags =	Net Int = 10969.5
----	---	--	---------	-------------------

Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

Ni	Conc = 0 mg/L Net Int %RSD = -76	Conc %RSD = Peak Int = 1566.00	Flags = b	Net Int = -98.7501
----	-------------------------------------	-----------------------------------	-----------	--------------------

Ni II 231.604 nm (mg/L)



Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	83753.9	4.000	3.972	-0.709
Standard 2	43082.7	2.000	2.023	1.14
Standard 3	23140.3	1.000	1.067	6.74
Standard 4	10969.5	0.5000	0.4842	-3.16
Blank	-98.7501	0	-0.04611	***

Calibration Coefficients	C0	C1	R
	-0.04137	4.791e-5	0.9992

Sample Blank [Sample Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

Ni	Conc = -0.04455 mg/L	Conc %RSD = -5.5	Flags = b	Net Int = -66.2425
	Net Int %RSD = -78	Peak Int = 1673.33		

1-1 Act Wt = 1.4546 Nom Wt = 25 Dil = 1

Ni	Conc = 19.90 mg/L	Conc %RSD = 2.57	Flags=b	Net Int = 24099.3
	Net Int %RSD = 2.54	Peak Int = 25831.3		

1-2 Act Wt = 1.6292 Nom Wt = 25 Dil =

Ni	Conc = 18.86 mg/L	Conc %RSD = 2.67	Flags = b	Net Int = 25617.4
	Net Int %RSD = 2.68	Peak Int = 27414.0		

2-1 Act Wt = 1.4015 Nom Wt = 25 Dil = 1

Ni	Conc = 11.68 mg/L	Conc %RSD = 1.86	Flags = b	Net Int = 13598.4
	Net Int %RSD = 1.87	Peak Int = 15321.3		

2-2 Act Wt = 1.5466 Nom Wt = 25 Dil = 1

Ni	Conc = 11.08 mg/L	Conc %RSD = 2.34	Flags =	Net Int = 14243.0
	Net Int %RSD = 2.36	Peak Int = 15978.0		

3-1 Act Wt = 1.4949 Nom Wt = 25 Dil = 1

Ni	Conc = 14.35 mg/L	Conc %RSD = 1.24	Flags = b	Net Int = 17846.5
	Net Int %RSD = 1.25	Peak Int = 19551.3		

3-3 Act Wt = 1.5796 Nom Wt = 25 Dil = 1

Ni	Conc = 14.03 mg/L	Conc %RSD = 2.46	Flags =	Net Int = 18429.1
	Net Int %RSD = 2.47	Peak Int = 20140.0		

Method Name :3

El	nm	Type	Mono	Order	Power (W)	Neb (L/min)	Height (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT Int. (V)	Int. (s)	Reps	Bkgd Corr	
Fe	II	238.204	An	1	2	1000	0.80	6.0	10.0	0.5	10	550	0.50	3	Fixed

Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Fe	Conc = 5.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 1035240
	Net Int %RSD = 2.37	Peak Int = 1048960		

Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Fe	Conc = 2.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 428150
	Net Int %RSD = 2.59	Peak Int = 438885		

Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Fe	Conc = 1.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 216089
	Net Int %RSD = 2.66	Peak Int = 225648		

Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Fe	Conc = 0.5000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 101647
	Net Int %RSD = 3.00	Peak Int = 110499		


Standard 5 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Fe	Conc = 0.1000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 25768.9
	Net Int %RSD = 1.59	Peak Int = 34231.3		

Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

Fe	Conc = 0 mg/L	Conc %RSD =	Flags = b	Net Int = -577.697
	Net Int %RSD = -12	Peak Int = 7847.33		

Fe II 238.204 nm (mg/L)



Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	1035240	5.000	4.979	-0.430
Standard 2	428150	2.000	2.047	2.35
Standard 3	216089	1.000	1.023	2.31
Standard 4	101647	0.5000	0.4705	-5.91
Standard 5	25768.9	0.1000	0.1041	4.06
Blank	-577.697	0	-0.02316	***

Calibration Coefficients	C0	C1	R
	-0.02037	4.829e-6	0.9997

Sample Blank [Sample Blank] Act Wt = 1 Nom Wt = 1 Dil = 1				
Fe	Conc = -0.003268 mg/L Net Int %RSD = 8.71	Conc %RSD = -46 Peak Int = 12056.7	Flags = b	Net Int = 3540.96
1-1 Act Wt = 0.8595 Nom Wt = 50 Dil = 2				
Fe	Conc = 38.29 mg/L Net Int %RSD = 1.47	Conc %RSD = 1.55 Peak Int = 80541.3	Flags =	Net Int = 71688.4
1-2 Act Wt = 0.7595 Nom Wt = 50 Dil = 2				
Fe	Conc = 36.28 mg/L Net Int %RSD = 2.75	Conc %RSD = 2.92 Peak Int = 69412.7	Flags =	Net Int = 60601.1
2-1 Act Wt = 1.4356 Nom Wt = 50 Dil = 2				
Fe	Conc = 9.059 mg/L Net Int %RSD = 1.87	Conc %RSD = 3.08 Peak Int = 15321.3	Flags = b	Net Int = 30475.0
2-2 Act Wt = 1.1685 Nom Wt = 50 Dil = 2				
Fe	Conc = 9.703 mg/L Net Int %RSD = 2.35	Conc %RSD = 2.70 Peak Int = 36548.0	Flags =	Net Int = 27021.2
3-1 Act Wt = 1.0860 Nom Wt = 50 Dil = 2				
Fe	Conc = 15.08 mg/L Net Int %RSD = 1.73	Conc %RSD = 1.91 Peak Int = 47198.7	Flags =	Net Int = 37449.4
3-2 Act Wt = 1.2581 Nom Wt = 50 Dil = 2				
Fe	Conc = 13.56 mg/L Net Int %RSD = 2.06	Conc %RSD = 2.37 Peak Int = 47578.7	Flags =	Net Int = 38865.5

Method Name :3


El	nm	Type	Mono	Order	Power (W)	Neb (L/min)	Height (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT (V)	Int. (s)	Reps	Bkgd Corr	
Ca	II	393.366	An	1	1	1000	0.80	6.0	10.0	0.5	10	400	0.50	3	Dyn

Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1				
Ca	Conc = 5.000 mg/L Net Int %RSD = 1.38	Conc %RSD = 0.00 Peak Int = 1250190	Flags =	Net Int = 1248720
Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1				
Ca	Conc = 2.000 mg/L Net Int %RSD = 2.00	Conc %RSD = 0.00 Peak Int = 535346	Flags =	Net Int = 534538
Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1				
Ca	Conc = 1.000 mg/L Net Int %RSD = 0.82	Conc %RSD = 0.00 Peak Int = 269316	Flags =	Net Int = 268688
Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1				
Ca	Conc = 0.5000 mg/L Net Int %RSD = 1.70	Conc %RSD = 0.00 Peak Int = 133854	Flags =	Net Int = 133321

Standard 5 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1				
Ca	Conc = 0.1000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 28148.8
	Net Int %RSD = 1.83	Peak Int = 28546.3		

Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1				
Ca	Conc = 0 mg/L	Conc %RSD =	Flags = b	Net Int = 2974.89
	Net Int %RSD = 0.92	Peak Int = 3378.33		

Ca II 393.366 nm (mg/L)



Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	1248720	5.000	4.955	-0.897
Standard 2	534538	2.000	2.095	4.74
Standard 3	268688	1.000	1.030	3.00
Standard 4	133321	0.5000	0.4878	-2.44
Standard 5	28148.8	0.1000	0.06658	-33.4
Blank	2974.89	0	-0.03424	***

Calibration Coefficients	C0	C1	R
	-0.04616	4.005e-6	0.9992

Sample Blank [Sample Blank] Act Wt = 1 Nom Wt = 1 Dil = 1				
Ca	Conc = 0.003954 mg/L	Conc %RSD = 27.6	Flags =	Net Int = 12511.8
	Net Int %RSD = 2.18	Peak Int = 12890.3		

1-1 Act Wt = 0.8595 Nom Wt = 50 Dil = 2				
Ca	Conc = 31.47 mg/L	Conc %RSD = 1.04	Flags =	Net Int = 80043.9
	Net Int %RSD = 0.88	Peak Int = 80537.0		

1-2 Act Wt = 0.7595 Nom Wt = 50 Dil = 2				
Ca	Conc = 29.89 mg/L	Conc %RSD = 2.45	Flags =	Net Int = 69184.6
	Net Int %RSD = 2.655	Peak Int = 69625.0		

2-1 Act Wt = 1.4356 Nom Wt = 50 Dil = 2				
Ca	Conc = 5.138 mg/L	Conc %RSD = 2.50	Flags = b	Net Int = 30929.7
	Net Int %RSD = 1.49	Peak Int = 31291.0		

2-2 Act Wt = 1.1685 Nom Wt = 50 Dil = 2				
Ca	Conc = 5.638 mg/L	Conc %RSD = 0.61	Flags =	Net Int = 28960.5
	Net Int %RSD = 0.35	Peak Int = 29361.7		

3-1 Act Wt = 1.0860 Nom Wt = 50 Dil = 2				
Ca	Conc = 7.179 mg/L	Conc %RSD = 2.40	Flags =	Net Int = 31978.5
	Net Int %RSD = 2.07	Peak Int = 32335.0		

3-2 Act Wt = 1.2581 Nom Wt =50Dil =2

Ca Conc = 7.041 mg/L Conc %RSD =1.57 Flags = Net Int = 34629.6
Net Int %RSD = 1.00 Peak Int = 35053.0

Method Name :3

El	nm	Type	Mono	Order	Power (W)	Neb (L/min)	Height (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT (V)	Int. (s)	Reps	Bkgd Corr
Mg II	279.553	An	1	2	1000	0.80	6.0	10.0	0.5	10	500	0.50	3	Dyn

Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Mg Conc = 1.000 mg/L Conc %RSD = 0.00 Flags = Net Int =1096130
Net Int %RSD = 2.15 Peak Int = 1102040

Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Mg Conc = 0.5000 mg/L Conc %RSD = 0.00 Flags = Net Int = 527422
Net Int %RSD = 1.57 Peak Int = 530994

Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Mg Conc = 0.1000 mg/L Conc %RSD = 0.00 Flags = Net Int = 112760
Net Int %RSD = 1.92 Peak Int = 114506


Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Mg Conc = 0.05000 mg/L Conc %RSD = 0.00 Flags = Net Int = 66841.7
Net Int %RSD = 0.75 Peak Int = 68448.0

Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

Mg Conc = 0 mg/L Conc %RSD = Flags = b Net Int =437.311
Net Int %RSD = 22.1 Peak Int = 1670.67

Mg II 279.553 nm (mg/L)




Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	1096130	1.000	1.007	0.719
Standard 2	527422	0.5000	0.4835	-3.30
Standard 3	112760	0.1000	0.1017	1.67
Standard 4	66841.7	0.05000	0.05939	18.8
Blank	437.311	0	-0.001759	***

Calibration Coefficients	C0	C1	R
	-0.00216	9.208e-7	0.9994

Sample Blank [Sample Blank] Act Wt = 1 Nom Wt = 1 Dil = 1															
Mg	Conc =0.5148 mg/L					Conc %RSD = 1.76					Flags =	Net Int =565843			
	Net Int %RSD = 1.75					Peak Int = 569502									
1-1 Act Wt =1.4546 Nom Wt = 25 Dil =1															
Mg	Conc =3.464 mg/L					Conc %RSD = 2.33					Flags=	Net Int = 781493			
	Net Int %RSD =045					Peak Int = 795689									
1-2 Act Wt = 1.6292 Nom Wt = 25 Dil = 1															
Mg	Conc = 3.697 mg/L					Conc %RSD = 3.30					Flags =	Net Int = 827479			
	Net Int %RSD =1.04					Peak Int = 832205									
2-1 Act Wt =1.4015 Nom Wt =25 Dil = 1															
Mg	Conc = 1.210 mg/L					Conc %RSD = 2.93					Flags = b	Net Int = 639519			
	Net Int %RSD = 2.38					Peak Int = 643404									
2-2 Act Wt = 1.5466 Nom Wt = 25 Dil = 1															
Mg	Conc = 1.769 mg/L					Conc %RSD = 2.41					Flags =	Net Int = 684710			
	Net Int %RSD = 2.44					Peak Int = 688857									
3-2 Act Wt = 1.3944 Nom Wt = 25 Dil = 1															
Mg	Conc =2.303 mg/L					Conc %RSD =2.25					Flags =	Net Int = 705331			
	Net Int %RSD = 2.47					Peak Int = 709511									
3-3 Act Wt = 1.5796 Nom Wt =25 Dil =1															
Mg	Conc = 2.653 mg/L					Conc %RSD =2.38					Flags =	Net Int = 747895			
	Net Int %RSD = 1.31					Peak Int = 752318									
Method Name :3															
El	nm	Type	Mono	Order	Power (W)	Neb (L/min)	Height (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT (V)	Int. (s)	Reps	Bkgd Corr	
Si	I	251.612	An	1	2	1000	0.80	6.0	10.0	0.5	10	500	0.50	3	Fixed
Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1															
Si	Conc = 4.000 mg/L					Conc %RSD = 0.00					Flags =	Net Int =40416.5			
	Net Int %RSD = 0.70					Peak Int = 42355.3									
Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1															
Si	Conc = 2.000 mg/L					Conc %RSD = 0.00					Flags =	Net Int = 20700.0			
	Net Int %RSD = 1.50					Peak Int = 22565.3									
Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1															
Si	Conc = 1.000 mg/L					Conc %RSD = 0.00					Flags =	Net Int = 11384.9			
	Net Int %RSD = 1.82					Peak Int = 13106.7									
Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1															
Si	Conc = 0.5000 mg/L					Conc %RSD = 0.00					Flags =0	Net Int = 5419.18			
	Net Int %RSD = 2.58					Peak Int = 7145.33									
Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1															

Si Conc = 0 mg/L Conc %RSD = Flags = b Net Int =291.363
 Net Int %RSD = 37.7 Peak Int = 2010.00

Si I 251.612 nm (mg/L)



Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	40416.5	4.000	3.979	-0.515
Standard 2	20700.0	2.000	2.006	0.308
Standard 3	11384.9	1.000	1.074	7.39
Standard 4	5419.18	0.5000	0.4769	-4.63
Blank	291.363	0	-0.03633	***

Calibration Coefficients	C0	C1	R
	-0.06549	1.001e-4	0.9992

Sample Blank [Sample Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

Si Conc =-0.01488 mg/L Conc %RSD =-19 Flags =b Net Int =517.454
 Net Int %RSD = 24.4 Peak Int = 2251.33

1-1 Act Wt =1.4546 Nom Wt = 25 Dil =1

Si Conc =6.173 mg/L Conc %RSD = 3.33 Flags= Net Int = 4103.45
 Net Int %RSD =3.12 Peak Int = 5752.00

1-2 Act Wt = 1.6292 Nom Wt = 25 Dil = 1

Si Conc = 5.706 mg/L Conc %RSD = 2.30 Flags = Net Int = 4229.84
 Net Int %RSD = 2.04 Peak Int = 5829.85

2-1 Act Wt =1.4015 Nom Wt =25 Dil = 1

Si Conc = 1.193 mg/L Conc %RSD = 3.93 Flags = Net Int = 1185.07
 Net Int %RSD = 3.38 Peak Int = 2830.00

2-2 Act Wt = 1.5466 Nom Wt = 25 Dil = 1

Si Conc = 1.254mg/L Conc %RSD = 4.41 Flags = Net Int = 1292.19
 Net Int %RSD = 4.44 Peak Int = 2965.33

3-1 Act Wt = 1.4949 Nom Wt = 25 Dil = 1

Si Conc =0.7304 mg/L Conc %RSD =3.15 Flags = Net Int = 953.515
 Net Int %RSD = 3.54 Peak Int = 2568.00

3-3 Act Wt = 1.5796 Nom Wt =25 Dil =1

Si Conc = 0.6514 mg/L Conc %RSD =2.38 Flags = Net Int = 928.381
 Net Int %RSD = 2.31 Peak Int = 2558.00

Method Name :W

El	nm	Type	Mono	Order	Power (W)	Neb (L/min)	Height (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT Int. (V)	Int. (s)	Reps	Bkgd Corr
As I	228.812	An	1	2	1000	0.80	6.0	10.0	0.5	10	700	0.50	3	Fixed

Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

As	Conc = 4.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int =65388.7
	Net Int %RSD = 1.39	Peak Int = 75174.7		

Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

As	Conc = 2.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 34146.0
	Net Int %RSD = 2.44	Peak Int =43705.3		

Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

As	Conc = 1.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 17525.6
	Net Int %RSD = 3.83	Peak Int = 26718.0		


Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

As	Conc = 0.5000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 7995.62
	Net Int %RSD = 4.58	Peak Int = 17240.7		

Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

As	Conc = 0 mg/L	Conc %RSD =	Flags = b	Net Int =233.629
	Net Int %RSD = 31.6	Peak Int = 9405.33		

As I 228.812 nm (mg/L)



Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	65388.7	4.000	3.966	-0.857
Standard 2	34146.0	2.000	2.056	2.78
Standard 3	17525.6	1.000	1.039	3.95
Standard 4	7995.62	0.5000	0.4569	-8.63
Blank	233.629	0	-0.01769	***

Calibration Coefficients	C0	C1	R
	-0.03197	6.114e-5	0.9992

Sample Blank [Sample Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

As	Conc =-0.02168 mg/L	Conc %RSD =-70	Flags =b	Net Int =168.349
	Net Int %RSD = 14.8	Peak Int = 9382.67		

1-1 Act Wt =1.4546 Nom Wt = 25 Dil =1

As	Conc =0.7537 mg/L	Conc %RSD = 10.7	Flags =	Net Int = 904.944
----	-------------------	------------------	---------	-------------------

Net Int %RSD =9.98		Peak Int = 9776.67	
1-2 Act Wt = 1.6292 Nom Wt = 25 Dil = 1			
As	Conc = 0.7990 mg/L	Conc %RSD = 5.30	Flags =
	Net Int %RSD =5.04	Peak Int = 10554.0	Net Int = 1039.22
2-1 Act Wt =1.4015 Nom Wt =25 Dil = 1			
As	Conc = 0.02635 mg/L	Conc %RSD = 10.3	Flags =
	Net Int %RSD = 9.84	Peak Int = 9418.67	Net Int = 212.334
2-2 Act Wt = 1.5466 Nom Wt = 25 Dil = 1			
As	Conc = 0.004426 mg/L	Conc %RSD = 9.41	Flags =
	Net Int %RSD = 8.44	Peak Int = 9399.03	Net Int = 192.667
3-1 Act Wt = 1.4949 Nom Wt = 25 Dil = 1			
As	Conc =0.06684 mg/L	Conc %RSD =8.15	Flags =
	Net Int %RSD = 7.54	Peak Int = 9460.52	Net Int = 253.515
3-3 Act Wt = 1.5796 Nom Wt =25 Dil =1			
As	Conc = 0.005447 mg/L	Conc %RSD =7.38	Flags =
	Net Int %RSD = 6.31	Peak Int = 9400.47	Net Int = 193.515

Method Name :W

El	nm	Type	Mono	Order	Power (W)	Neb (L/min)	Height (mm)	Plasma (L/min)	Aux (L/min)	Pump (RPM)	PMT Int. (V)	Reps	Bkgd Corr
Cu I	324.754	An	1	1	1000	0.50	6.0	10.0	0.5	10	500	0.50	3 Dyn

Standard 1 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Cu	Conc = 3.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int =54379.0
	Net Int %RSD = 1.44	Peak Int = 58069.3		

Standard 2 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Cu	Conc = 2.000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 36193.9
	Net Int %RSD = 1.23	Peak Int =39798.7		

Standard 3 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Cu	Conc = 1.500 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 27294.0
	Net Int %RSD =2.97	Peak Int = 30877.3		


Standard 4 [Cal. Std] Act Wt = 1 Nom Wt = 1 Dil = 1

Cu	Conc = 0.5000 mg/L	Conc %RSD = 0.00	Flags =	Net Int = 18835.3
	Net Int %RSD = 0.69	Peak Int = 12570.0		

Blank [Cal. Blank] Act Wt = 1 Nom Wt = 1 Dil = 1

Cu	Conc = 0 mg/L	Conc %RSD =	Flags = b	Net Int = 679.034
	Net Int %RSD = 14.9	Peak Int = 4350.67		

Cu I 324.754 nm (mg/L)



Standard	Counts	Real Conc	Calc Conc	% Diff
Standard 1	54379.0	3.000	3.005	0.154
Standard 2	36193.9	2.000	1.991	-0.458
Standard 3	27294.0	1.500	1.495	-0.355
Standard 4	18835.3	1.000	1.023	2.31
Standard 5	9017.26	0.5000	0.4758	-4.84
Blank	679.034	0	0.01094	***

Calibration Coefficients	C0	C1	R
	-0.02692	5.575e-5	0.9998

Sample Blank [Sample Blank]	Act Wt = 1	Nom Wt = 1	Dil = 1		
Cu	Conc = -0.003589 mg/L	Conc %RSD = -68	Flags = b	Net Int = 418.503	
	Net Int %RSD = 10.4	Peak Int = 4046.67			
1-1 Act Wt = 0.8594 Nom Wt = 50 Dil = 1					
Cu	Conc = 0.3255 mg/L	Conc %RSD = 6.78	Flags =	Net Int = 524.449	
	Net Int %RSD = 5.98	Peak Int = 4328.67			
1-2 Act Wt = 0.7595 Nom Wt = 50 Dil = 1					
Cu	Conc = 0.3626 mg/L	Conc %RSD = 7.30	Flags =	Net Int = 522.898	
	Net Int %RSD = 6.04	Peak Int = 4216.67			
2-1 Act Wt = 1.4356 Nom Wt = 50 Dil = 1					
Cu	Conc = 0.3343 mg/L	Conc %RSD = 6.20	Flags =	Net Int = 596.263	
	Net Int %RSD = 5.84	Peak Int = 4302.00			
2-2 Act Wt = 1.1685 Nom Wt = 50 Dil = 1					
Cu	Conc = 0.3621 mg/L	Conc %RSD = 5.41	Flags =	Net Int = 570.279	
	Net Int %RSD = 5.44	Peak Int = 4361.33			
3-1 Act Wt = 1.086 Nom Wt = 50 Dil = 1					
Cu	Conc = 0.6141 mg/L	Conc %RSD = 6.15	Flags =	Net Int = 663.326	
	Net Int %RSD = 6.54	Peak Int = 4313.33			
3-2 Act Wt = 1.2581 Nom Wt = 50 Dil = 1					
Cu	Conc = 0.6831 mg/L	Conc %RSD = 5.38	Flags =	Net Int = 726.802	
	Net Int %RSD = 5.31	Peak Int = 4320.67			

样品检测结果

单位：mg/kg

元素	样品 1	样品 2	样品 3
Cu	0.34	0.35	0.65
Pb	未检出	未检出	未检出
Fe	37.29	9.38	14.32
Ni	19.38	11.38	14.19
Na	51.68	21.15	20.90
Ca	30.68	5.39	7.11
Mg	3.58	1.49	2.48
Al	147.15	4.13	12.37
As	0.78	未检出	未检出
Si	5.94	1.22	0.69
V	0.66	1.86	5.43