

Analysis of frankincense essential oil compositions by GC-MS

1 Instruments and reagents

GC-MS3100 gas chromatograph/mass spectrometer (EAST & WEST ANALYTICAL INSTRUMENTS, INC); Equity-5 (30m×0.25mm×0.25um) silica capillary column; 1 µL micro injector;

Sample: frankincense essential oil

2 Analysis conditions

MS conditions: EI source; ion source temperature: 200 °C; electron energy: 70eV; scanning mode: full scan; scanning range: 28.5u~400u; scanning period: 0.6 s; interface temperature: 250 °C; multiplier high voltage: 1200 V.

GC conditions: Equity-5 (30m×0.25mm×0.25um) silica capillary column; injection port temperature: 250 °C; split sampling; sample volume: 0.08 µL; column flow: 1 mL/min; split ratio: 70:1; precolumn pressure: 50 kPa; purge rate: 2ml/min; column temperature program: hold at 50 °C for 1min, ramp up to 250 °C at 5 °C/min and then hold for 10 min.




Figure 1 TIC chromatogram of frankincense essential oil




Figure 2 Partially magnified chromatogram

Table 1 Analysis results of frankincense essential oil composition

Peak No.	Retention time, min	Name	CAS No	Molecular formula	relative amount, %	similarity degree, %
1	7.15	1-Isopropyl-4-methylbicyclo[3.1.0]hex-2-ene	58037-87-9	C ₁₀ H ₁₆	64.41	95
2	7.30	Bicyclo[3.1.1]hept-2-ene, 2,6,6-trimethyl-, (1R)-	7785-70-8	C ₁₀ H ₁₆	5.43	93
3	7.55	Bicyclo[3.1.0]hex-2-ene, 4-methylene-1-(1-methylethyl)- (2,4(10)-Thujadiene)	36262-09-6	C ₁₀ H ₁₄	0.21	83
4	7.68	4-Carene, (1S,3R,6R)-(-)-	5208-49-1	C ₁₀ H ₁₆	0.06	74
5	8.40	Phellandrene	555-10-2	C ₁₀ H ₁₆	4.92	90
6	8.50	beta-Pinene	127-91-3	C ₁₀ H ₁₆	0.29	83
7	8.86	1,6-Octadiene, 7-methyl-3-methylene- (Myrcene)	123-35-3	C ₁₀ H ₁₆	0.67	
8	9.30	1,3-Cyclohexadiene, 2-methyl-5-(1-methylethyl)-	99-83-2	C ₁₀ H ₁₆	2.14	91
9	9.49	3-Carene	13466-78-9	C ₁₀ H ₁₆	4.14	
10	9.66	Cyclohexene, 1-methyl-4-(1-methylethylidene)-	586-62-9	C ₁₀ H ₁₆	0.21	85
11	9.91	Benzene, 1-methyl-4-(1-methylethyl)- (p-Cymene)	99-87-6	C ₁₀ H ₁₄	2.40	90
12	10.05	D-Limonene	5989-27-5	C ₁₀ H ₁₆	2.11	88
13	10.28	1,3,6-Octatriene, 3,7-dimethyl-, (E)-	3779-61-1	C ₁₀ H ₁₆	0.43	91
14	10.60	1,3,6-Octatriene, 3,7-dimethyl-, (Z)-	3338-55-4	C ₁₀ H ₁₆	0.17	81

Peak No.	Retention time, min	Name	CAS No	Molecular formula	relative amount, %	similarity degree, %
15	10.95	1,4-Cyclohexadiene, 1-methyl-4-(1-methylethyl)-	99-85-4	C ₁₀ H ₁₆	0.47	85
16	11.26	5-Isopropyl-2-methylbicyclo[3.1.0]hexan-2-ol cis-Sabinenehydrate	15826-82-1	C ₁₀ H ₁₈ O	0.12	82
17	11.87	2-Carene	554-61-0	C ₁₀ H ₁₆	0.18	86
18	12.22	Terpineol, cis-	7299-41-4	C ₁₀ H ₁₈ O	0.20	80
19	12.78	Thujone	546-80-5	C ₁₀ H ₁₆ O	0.16	86
20	14.67	1-Methyl-4-isopropyl-1-cyclohexen-4-ol (4- Terpineol)	562-74-3	C ₁₀ H ₁₈ O	0.76	84
21	15.28	Anisole, p-allyl-	140-67-0	C ₁₀ H ₁₂ O	2.61	92
22	20.55	Tricyclo[4.4.0.0(2,7)]dec-3-ene, 8-isopropyl-1,3-dimethyl-	3856-25-5	C ₁₅ H ₂₄	0.10	74
23	20.82	Cyclobuta[1,2:3,4]dicyclopentene, decahydro-3a-methyl-6-methylene-1-(1-methylethyl)-,	5208-59-3	C ₁₅ H ₂₄	0.46	88
24	21.19	Benzene, 1,2-dimethoxy-4-(2-propenyl)-	93-15-2	C ₁₁ H ₁₄ O ₂	0.26	84
25	21.26	Humulen-(v1)	N/A	C ₁₅ H ₂₄	0.03	74
26	22.04	Germacrene D	23986-74-5	C ₁₅ H ₂₄	0.04	76
27	22.12	Bicyclo[3.1.1]hept-2-ene, 2,6-dimethyl-6-(4-methyl-3-pentenyl)-	17699-05-7	C ₁₅ H ₂₄	0.04	78
28	23.43	Naphthalene, 1,2,4a,5,6,8a-hexahydro-4,7-dimethyl-1-(1-methylethyl)-	483-75-0	C ₁₅ H ₂₄	0.04	80
29	23.50	Isocaryophyllene	N/A	C ₁₅ H ₂₄	0.05	81
30	24.48	Cadina-1(10),4-diene	483-76-1	C ₁₅ H ₂₄	0.09	70
31	24.68	VERIDIFLOROL	N/A	C ₁₅ H ₂₆ O	0.55	68
32	25.17	Benzene, 1,2,3-trimethoxy-5-(2-propenyl)-	487-11-6	C ₁₂ H ₁₆ O ₃	0.07	81
33	31.03	Bicyclo[3.1.0]hex-2-ene, 2-methyl-5-(1-methylethyl)- (3- Thujene)	2867-05-2	C ₁₀ H ₁₆	0.05	76
34	31.50	Isopropyl Myristate	110-27-0	C ₁₇ H ₃₄ O ₂	6.14	88